



SAMARIT

safe patient moving & handling

SAMARIT[™]
Rollbord[™]
ECOLITE

For
medium
use!

SAMARIT "no lift" Transfer Systems

Rollbord[™] **ECOLITE**

The safe way to transfer patients

A legal requirement in some countries, taken for granted in others – working conditions should be as comfortable as possible. Healthcare workers should avoid lifting when transferring patients. **SAMARIT Products** assist with transfer, without the need to lift. The products are designed to facilitate the transfer without using a lot of effort, thus making the work easier for the carers and kinder to their backs.

For medium use!

3 inexpensive EcoLite models for medium use.

EcoLite 180

foldable

For radiology, operating theatre, intensive care units and all special wards.



*SWL: 160 kg / 353 lbs

Weight: 2.8 kg
Size: 180 x 50 cm, 89 x 50 cm
70 x 20 inch / 35 x 20 inch

Order No. 440.1180
UDI/GTIN: 7640144460327

EcoLite 120

foldable

For transfer of injured patients.



*SWL: 150 kg / 331 lbs

Weight: 1.4 kg
Size: 130 x 40 cm, 65 x 40 cm
51 x 16 inch, 25 x 16 inch

Order No. 440.1120
UDI/GTIN: 7640144460334

EcoLite 90

non-foldable

For urology operating theatres, gynaecology operating theatres, wake-up stations and emergencies.



*SWL: 150 kg / 331 lbs

Weight: 1.3 kg
Size: 90 x 50 cm, 35 x 20 inch

Order No. 440.1090
UDI/GTIN: 7640144460341

Environmentally friendly material SAMARIT** transfer aids are made of ecologically safe materials. They can be disposed of in an environmentally friendly manner and if a fire should occur they do not produce toxic fumes.

** Silicone sliding fabric around a polyethylene foam core, coated polyester cloth, PE ends.

* SWL = Safe Weight Load

EU-SRN: CH-MF-000016280
CHRN: CHRN-MF-20000184

www.samarit.com

Easy cleaning

The Rollbord can be cleaned best with *SamaClean Wipes* or any hospital disinfectant. For anticipated soiled conditions use the Rollbord disposable cover. It covers ends of Rollbord and allows any residue to be folded untouched into the cover, leaving a clean surface. Repairs can only be carried out by authorized personnel (warranty).

High quality standards

SAMARIT transfer aids have:

- CE mark
- Patents and international design patents
- 2 years warranty



SAMARIT Rollbord Hightec



SAMARIT Transglide



SAMARIT Rollbody



SAMARIT Nursing Aids

Download your language from our website or ask your local distributor.



SAMARIT Medical AG

Gewerbestr. 12

CH-8132 Egg, Schweiz

Tel.: +41 44 918 10 11 Fax: +41 44 918 24 39

e-mail: info@samarit.com - www.samarit.com

